



PRODUCT CATALOG FOR LAB



Dental 
Save Money. Business Better!



Guangzhou YiLinK Medical Equipment Technology Co.,Ltd
Tel: 020-32208271
Web: WWW.yilink-medical.com
Fax: 020-32208069
Add: Building F 405, Software Science Park, No.11 Caipin Road, Luogang Science City, Huangpu District, Guangzhou.





AESTHETIC AND REALISTIC resembling natural teeth PMMA

- ✓ MONOCOLOR
- ✓ MULTILAYER
- ✓ FLEXIBLE
- ✓ MULTI-MIX COLOR
- ✓ CLEAR
- ✓ PINK



ONE PRICE FOR ALL SIZES
Save Money. Business Better!




Monocolor PMMA

Yilink Monocolor PMMA is used for temporary restorations, 16 colors available, a perfect fitting temporary is quickly produced. Various available geometries and sizes makes it compatible with all kinds of milling machines. Yilink Temp can also be used as a fit and function check before the final restoration is produced. We recommend this PMMA for machine calibration as well.

Recommended Indications



Colors



Thickness

10mm 12mm 14mm 16mm 18mm 20mm 22mm 25mm

Material Characteristic:

Thermal Expansivity (25-500 c)[10/K]	10.5±0.5
Bending Strength [MPa]	400±60
Chemical Solubility [ug/cm ²]	≤ 40
Density[g/cm ³]	2.47±0.05

Advantage

1. Various colors and geometries
2. Long-term color stability and easy to polish
3. Biocompatible because of low residual monomer



Material Characteristic:

Flexural Strength	≥ 120[MPa]
Density	≥ 1.18 [g/m ³]
Vickers hardness	≥ 20[HV1]
Water absorption	≤ 25[μ g/mm ³]
Solubility	≤ 5[μ g/mm ³]
Residual monomer	≤ 0.5[%]

Multilayer PMMA

The aesthetic solution for temporary restorations

Yilink Multilayer PMMA is the aesthetic solution for temporary restorations . Thanks to the smooth layertransition and natural color gradient , only polishing is necessary to finalize long term temporary crown andbridges . Available in 16 shade as well as various geometries and sizes ,YilinkTemp iscompatible with all kinds of milling machines .

Recommended Indications



Colors



Thickness

10mm12mm14mm16mm18mm20mm22mm25mm

Advantage

- 1.Esthetic,natural results only through polishing
- 2.Various colors and geometries
- 3.High strength and biocompatible for long term temporaries



Material Characteristic:

Flexural Strength	≥ 100[MPa]
Density	≥ 1.18 [g/m ³]
Vickers hardness	≥ 20[HV1]
Water absorption	≤ 40[μ g/mm ³]
Solubility	≤ 5[μ g/mm ³]
Residual monomer	≤ 0.5[%]

Multi-Mix Color PMMA


Multi-mix color full denture dental pmma

Yilink multi mix color PMMA is used for full denture. In terms of color, it can provide 16 colors, and can match a variety of sizes to meet the needs of different mainstream cutting systems

Recommended Indications



Colors



Thickness

30mm 40mm

Advantage

- 1.Stable, denture won't fall apart from gingival base
- 2.Save time and labor cost
- 3.Excellent long-term stability and esthetics